




OTS ELITE

AUTO-TRACKING RADIOGRAPHIC SYSTEM

Excellence in Workmanship, Best-in-Class of High Performance Medical Systems

The OTS Elite delivers maximum capability and is designed for use in hospitals high-volume orthopedics, and free-standing ERs. Loaded with convenience features like rotational receptor trays, patient overhead grip and table-top lock release and elevation handle, the OTS Elite stands out in the class of high performance medical systems by offering the only 5 year parts warranty of its kind! Technologist-friendly features, such as color-coded annotations and the ability to turn on or off the auto-tracking with the touch of a button make positioning fast, improving patient throughput. Improve your patients' experience... choose the AMRAD[®] Medical OTS Elite.





OTS

Optimally balanced ceiling mounted tube stand makes for featherweight, easy movement

Ceiling mounted tube stand for maximum flexibility

Recommended ceiling height 9' to eliminate any positioning challenges

Perfect for facilities who want improved patient throughput

Auto-tracking to elevating table and wallstand

J1000 WALLSTAND

Ergonomic lock-release for ease of positioning

PA Hand grips for patient stability

Overhead Patient Grip for positioning ease

Conveniently placed chin rest ensures patient comfort and aligns the body to the receptor for upright exams



Rotational receptor trays help to protect your DR investment



PL

TE

S222 TABLE

650 lb. patient weight capacity

Whisper quiet for patient comfort when they need it the most

Programmable working height to fit any situation


Recessed four pedal design provides easy operation and patient safety

Smart Investment

The OTS Elite offers equal features and earns equal reimbursement to all

of more costly competition and all of this at 5 times the warranty.

Smart Decisions equal smart investments. Be Smart!





www.amradmedical.com

SUMMIT INDUSTRIES, LLC 7555 N. Caldwell Avenue, Niles, Illinois 60714 (p) 773-353-4000 | (p) 1-800-729-9729

