

# PERFORMANCE RUNS IN THE FAMILY.

A Choice of Models Delivers  
the Full Power of Digital.



# POWERFUL. DIGITAL. RELIABLE.

With the DRX Detector, Carestream introduced the world's first wireless, cassette-sized solution. Today, our family of DRX Core and DRX Plus Detectors carry on this tradition of innovation.

Our affordable DRX Core Detector converts the X-ray equipment you're using right now, upgrading it to the benefits of full digital.

The DRX Plus, with X-Factor shareability, is thinner and lighter for easier handling. And, it offers an IPX 57 water-resistance rating for greater reliability.

Finally, the DRX 2530C delivers X-Factor shareability in a smaller format – easy to position in incubator X-ray trays.

# GET READY TO MEET THE FAMILY.

# THE INDUSTRY BENCHMARK FOR DIGITAL IMAGE CAPTURE.

Carestream's DRX Detectors let you upgrade from film or CR to full-digital radiography\* – quickly, easily and affordably. These solutions provide high-quality X-ray images in seconds, for rapid, confident diagnoses and earlier commencement of treatment. Plus, we offer a wide array of Detector models, so you can select precisely the right ones for your facility.

DRX is a proven technology – small wonder there are over 12,000 of our Detectors in use worldwide.

\* As part of an integrated system including an optional Carestream console, preloaded software and charger.

## A Family of Options.

DRX Detectors are available with:

- Gadolinium scintillator for general radiology imaging.
- Cesium Iodide scintillator for low-dose applications.
- 35x43cm, 43x43cm and 25x30cm sizes.
- Wireless or fixed design.



### Carestream Radiography Systems.

Our complete line of innovative DR image-capture solutions, designed around the DRX Detector, includes full radiology rooms and mobile x-ray systems.



## Digital Made Affordable: The DRX Core Detector

The DRX Core lets you say goodbye to film or CR, gaining the benefits only DR can provide, without breaking your budget. The DRX Core includes:

- The ability to upgrade your existing systems without the high cost of equipment replacement.
- The superb image quality of full digital.
- Wireless connectivity (optional) to eliminate the hassles and hazards of cables and reduce the risk of infection.
- IPX 44 water-resistance rating, for enhanced reliability.
- Beam-detect mode, for cable-free, automatic capture activation.
- Carbon-fiber twill panel for improved handling.
- Detector-status LED indicators.



## Top-of-the-Line Performance: The DRX Plus Detector

Our DRX Plus delivers the power of X-Factor thinking in a thinner, lighter design for easier handling. In addition to all DRX Core features, the DRX Plus also offers:

- The ability to share detectors between any DR room or DR mobile system.
- Standard wireless connectivity to eliminate the hassles and hazards of cables and reduce the risk of infection.
- The power to upgrade virtually any manufacturer's system to DR
- IPX 57 water-resistance rating for superb reliability – submersible in up to a meter of water without failure.
- Battery hot-swap and a separate charger, allowing the detector to be used while a spare battery is being charged.



## Powerful Performance in a Smaller Package: The DRX 2530C Detector

The DRX 2530C delivers the power of X-Factor thinking in a smaller-format design, delivering:

- Fast, easy positioning in pediatric incubator trays.
- A Cesium iodide (CsI) scintillator for dose-sensitive pediatric applications.
- Compatibility with the DRX-Revolution small-format detector holder.
- Faster orthopaedic tabletop imaging due to its easy handling.
- An ergonomically designed handle for fast and safe positioning.
- Full compatibility with the same DRX battery charger.

DRX-Ascend



DRX-Revolution



DRX-Mobile  
Retrofit Kit



DRX-Transportable/  
Field Portable



NO MATTER WHICH ROAD YOU'RE ON,  
WE'LL GET YOU WHERE  
YOU NEED TO GO.

RIGHT FOR TODAY.  
READY FOR TOMORROW.



### Right for Today. Ready for Tomorrow.

Carestream is ready to help you plan the most effective route to your X-ray imaging future. Our scalable equipment design and modular components mean high performance today, along with easy, affordable upgrades for years to come.

Let's plan your equipment migration together – you'll gain the confidence that your current technology investment will continue to pay dividends well into your future.

### Protect Against Drops and Other Damage.

All DRX Detectors are eligible for coverage under the DRX Accident Protection Plan. It can be added to your standard Carestream DRX Warranty or Service Agreement, providing replacement of your detector at a fraction of the full cost, should damage occur during normal clinical use. For you, this means peace-of-mind, one-call resolution and fast replacement.

### A Community of Service and Support

For dependable service, look to our Customer Success Network.

We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.



[carestream.com](http://carestream.com)

