

# KMC 950

## High Performance Surgical C-arm

The KMC-950 is revolutionary in its versatility, efficiency and exceptional image resolution. Our improved touch screen control panel boasts expanded memory, providing an abundance of customizable menus, unmatched control versatility and virtually maintenance-free operation.



## Specifications



12kW, High Frequency Generator



135° Wide Orbital Rotation with large SID distance



Easy Archiving and Documentation

- DICOM 3.0
- CD/DVD Burner
- DICOM Paper Print Support
- USB Port



1Kx1k High Resolution CCD Camera



Touch Screen C-Arm control



High Performance Rotating Anode X-Ray Tube



High Resolution B/W 19" LCD Medical Grade Monitors

## Smart CX-View



Database Management

- More convenient "search&save" capabilities
- Manage patient information
- Before & after surgery comparison



Versatile Diagnostic Tools

- Zoom, Contrast, Reverse, Invert
- Filtering, Edge Enhancement
- Annotation



Unlimited Image Memory

- Unlimited patient capacity
- Unlimited study images



PACS connectivity with DICOM 3.0

- Store and Print, Modality Worklist Management
- MPPS, Storage Commitment, and Import & Export



## Option

- DSA (Digital Subtraction Angiography) with roadmap
- DAP (Dose Area Product)
- Virtual colimator
- Premium Cart(Folding and elevating twin 19" LCD monitor)
- Thermal printer



**GEMSS MEDICAL SYSTEMS CO., LTD.**

61, Dunchon-daero 541 beon-gil, Jungwon-gu,  
Seongnam-si, Gyeonggi-do, Korea  
Phone: +82-31-764-7321  
Fax: +82-31-764-7324  
[www.gemss-medical.com](http://www.gemss-medical.com)

The Power of Passion

**GEMSS**