


# EXOTIC DIGITAL X-RAY



Preview monitor is an upgrade from SimonDR

**SimonDR**

# EXOTIC DIGITAL X-RAY

A premium x-ray system designed by Quantum Medical Imaging combined with digital x-ray integration by SimonDR, Inc. brings a no compromise digital radiography system that is easy to use. The SimonDR DRV System is powered by WhiteCap Image viewing and processing software. It is designed to be simple and easy to use.

## Exotic Features Include:


- \* SimonDR DRV high-resolution CCD based direct digital x-ray system.
- \* 5.0 line pairs per millimeter.
- \* 20 million pixel sensor.
- \* FOV 12x17
- \* CsI (TI) Scintillator
- \* Integrated WhiteCap software generator control.
- \* Made and supported in the U.S.A.
- \* 32kw - 25ma - 500ma high-frequency standard power.
- \* Higher ma values mean faster exposure times.
- \* 200,000 heat-unit tube with a standard fractional focal spot (.6mm / 1.5mm).
- \* 60 inch four-way floating table top with 30 inches of longitudinal float to accommodate larger animals.
- \* Included "Uro Track" for easy liquid containment.
- \* Patient tie-downs.
- \* Recessed foot-control operated floating table-top with premium fail-safe brakes.
- \* 650 pound patient capacity tabletop.
- \* Laser guided collimator for ultimate accuracy.
- \* Pre-set anatomical exposure factors capable.

## Standard Warranty:

2 year receptor and standard 5 year limited x-ray generation warranty.

## Options Include:

- \* Stored energy gives you the ability to use high-frequency where incoming power is limited.
- \* Shorter float for very small areas.
- \* Tubestand-mounted preview monitor.



800.835.3852  
www.simondr.com

**Simon***dr*