



CAPITAL X-RAY, INC.

Carestream PRO

--

DR Retrofit

Carestream  
Branded  
DR Retrofit



Carestream

**Carestream PRO DR Retrofit System**

This state-of-the-art system provides the ideal balance of price & performance for a wide variety of imaging facilities -- small hospitals, orthopedic practices, imaging centers, urgent care facilities -- the list goes on and on!

Choose between a 35x43cm (14x17") wireless or tethered panel, OR a 43x43cm (17x17") fixed panel for a fast & easy cost-effective upgrade to digital. It's as easy as sliding the detector into your existing radiographic equipment!

**Increase Workflow While Enhancing Patient Care**

- Beam-detect mode allows for no cable connection to the generator -- the image capture begins automatically when the detector senses the x-ray exposure.
- Wireless option eliminates the hassle of cables and minimizes interference with ICU equipment.
- Available in Gadolinium or Cesium allowing you to further monitor and reduce the radiation dose to patients.

**Includes the Following:**

- Your choice of detector (Gadox or Cesium)
- Your choice of console (desktop CPU or laptop configuration)
- Pre-loaded with Carestream VXvue software
- ImageSuite Mini-PACS module
- Three-Year DR manufacturer's warranty



Contact us today  
to discuss  
the best  
configuration for your  
facility!

1-800-239-9729

Promotional Pricing  
Available

For more information, contact us today!

T: 800-239-9729 [www.capitalxray.com](http://www.capitalxray.com)

Image  
Is  
Everything.